

ALL - BOND UNIVERSAL



« Darf ich mich Ihnen vorstellen:



Ich bin hydrophil bei der Applikation für eine ideale Oberflächenbenetzung und hydrophob nach der Polymerisation für einen optimalen Verschluss und optimale Resistenz.

Ich brauche keinen Aktivator oder Primer, d.h. Sie sparen Zeit und Geld.

Ich kann sowohl nach dem Self-Etch, Selective Etch oder Total Etch Konzept eingesetzt werden.

Ich enthalte MDP & BPDM, ideal für Bonding auf Zirkonoxid, Aluminiumoxid und Metallen.

Ich habe eine Filmstärke von unter 10 Mikron, habe also keine negativen Auswirkungen auf die Passung von Rekonstruktionen.

Ich bin kompatibel mit allen dual- oder selbsthärtenden Zementen.

Ich bin der ideale Primer für eine optimale Adhäsion aller (in)direkten Rekonstruktionen.


Ich bin die erste Wahl für jede Stuhlassistenz, denn von jetzt ab bin nur noch ich für alle Bondingverfahren erforderlich. »

All-Bond Universal™ ist das einzige universelle lichthärtende Einflaschen-Adhäsivsystem, das ohne den Einsatz von Aktivatoren für alle gängigen dualhärtenden Zemente geeignet ist.

Bringing **Science** to
the *Art* of Dentistry™

« ... Und selbstverständlich sind meine Aussagen wissenschaftlich untermauert »

Microtensile Bond Strength on Dentin**




Dr. A. Wagner, Dr. R. Belli, Dr. J. Zorzin,
Prof. A. Petschelt, Prof. U. Lohbauer
(University of Erlangen, Germany)

Academy of Dental Materials 2013 Conference – Dent Mater 2013; 29s:e87-88




Dentin Bond Strength Comparison: Indirect Restorations*



ALL-BOND UNIVERSAL was cemented with DUO-LINK™, and ScotchBond™ Universal was cemented with RelyX™ ARC. Cements were self-cured.

*data on file



DUO-LINK was used with ALL-BOND UNIVERSAL, and Multilink® Automix was used with Multilink Primer A and B. Relyx Unicem 2 is a self-adhesive resin cement. Cements were self-cured. DENTIN WAS NOT ETCHED.

Trademarks are property of their respective owners.

BESTELLINFORMATION

BIB-7202P All-Bond Universal 6ml Flasche

CHF 153.45

Alle Preise zzgl. MwSt. und Versandkosten. Preisänderungen vorbehalten.

ms-dental ⁺

MS Dental AG · Dentalprodukte · Produits dentaires · Juraweg 5 · CH-3292 Busswil
T +41 32 387 38 68 · F +41 32 387 38 78 · info@msdental.ch · www.msdental.ch